

 <b>Fraunhofer</b> MEVIS	GAMMASTAR Server Release Letter	V1	Document-ID GSTAR_SERVER Last Change on 2025-09-11	Last Change by dhoi Page 1 / 3
--	------------------------------------	----	---	---

Last Change	Review	Approval

## Document History

Version	Date	Changes	JIRA-Ticket	Name
V1	2025-09-11	Release of v1.0 of gammaSTAR Server	n/a	dhoi

 <b>Fraunhofer</b> MEVIS	GAMMASTAR Server	Document-ID GSTAR_SERVER	Last Change by dhoi
	Release Letter	V1	Last Change on 2025-09-11

## 1 Introduction and Overview

This document summarizes the delivery of the **gammaSTAR Sequence Server** for receiving raw MR sequence events and forward it to a simulation tool.

Contact person at Fraunhofer MEVIS: Daniel Christopher Hoinkiss (Daniel.hoinkiss@mevis.fraunhofer.de)

## 2 Release Summary

<b>Software Component</b>	gammaSTAR Server
<b>Version</b>	<a href="https://gitlab.fme.lan/gammastar/gs-py_backend/-/tree/gammastar_release/sequenceServer?ref_type=heads">https://gitlab.fme.lan/gammastar/gs-py_backend/-/tree/gammastar_release/sequenceServer?ref_type=heads</a> <ul style="list-style-type: none"> <li>SHA: ae991d132b189b043cc9d2717a551288e9745af1</li> </ul>
<b>Embedding Product(s)</b>	Not applicable
<b>Purpose of this Release</b>	Provision of the gammaSTAR Server v1.0 software under AGPLv3-clause license for research and educational purposes. Not for clinical use.
<b>OS and Compiler Details</b>	OS: Windows 11 Python Version: 3.12
<b>Embedded Libraries</b>	<p>The following libraries are used by the software. Detailed license information is contained in the delivered "third_party_licenses" folder and in the NOTICE file.</p> <ul style="list-style-type: none"> <li>Websockets - BSD License</li> <li>numpy - BSD License</li> <li>matplotlib - Matplotlib License (BSD-style)</li> <li>isimrmrd - ISMRMRD SOFTWARE LICENSE JULY 2013</li> <li>h5py - BSD License</li> <li>pydicom - MIT License</li> <li>fastapi - MIT License</li> <li>MRzero - AGPL-3.0 license</li> </ul>
<b>Functionality</b>	The software acts as a fastAPI API that can receive gammaSTAR raw representation and forwards them to MR Zero. It forwards the incoming simulated or measured raw data to the gammaSTAR image reconstruction.
<b>Regulatory Compliance</b>	<p>THIS SOFTWARE IS PROVIDED "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED.</p> <p>In particular, Fraunhofer MEVIS does not ensure compliance to medical product regulations for the release.</p>
<b>Safety Classification</b>	None

 <b>Fraunhofer</b> MEVIS	GAMMASTAR Server Release Letter	Document-ID GSTAR_SERVER	Last Change by dhoi
	V1	Last Change on 2025-09-11	Page 3 / 3

<b>Risk Measures</b>	No risk measures defined.
<b>Known Nonconformities and Restrictions</b>	The software is not qualified for use as a medical product or as part thereof. No bugs or restrictions are known.